

B.O.D.seed for Accurate B.O.D. Testing



Description

B.O.D.seed increases accuracy and consistency in B.O.D. testing and eliminates costs associated with collecting and maintaining an acclimatized seed.

Each **B.O.D.seed** capsule contains specialized microbial cultures to provide a broad range of organisms suitable for most types of industrial and municipal waste water.

One capsule of **B.O.D.seed** provides enough acclimatized seed for up to 250 B.O.D. tests daily, at a cost of pennies per test and in strict accordance with "Standards Methods" as established by the U.S.E.P.A.

Specifications

Form: Free-flowing granular powder in a capsule, 50 capsules in a bottle.
Color: Brown, capsule is clear

B.O.D.seed for Accurate B.O.D. Testing

**ISO 9001: 2000 Certified
Quality System**



**MANUFACTURED IN USA SINCE
1987**

CERTIFIED PATHOGEN FREE

